


YASARGIL MINI VESSEL CLIP

permanent; 5mm; opening 3.5mm

130 gms; straight



SKU: 244-473-05

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE716
AUTHORS	YASARGIL
Medicon	58.57.40
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE716
AUTHORS	YASARGIL
Medicon	58.57.40

SPECIFICATIONS

[permanent](#), [straight](#)