


YASARGIL MINI VESSEL CLIP

permanent; 4.5mm; opening 6mm

130 gms; angled



SKU: 244-477-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE806
AUTHORS	YASARGIL
Medicon	58.57.48
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE806
AUTHORS	YASARGIL
Medicon	58.57.48

SPECIFICATIONS

[angled](#), [permanent](#)