


YASARGIL MINI VESSEL CLIP

temporary; 3mm; opening 3.3mm

90 gms; straight



SKU: 244-560-03

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE681
AUTHORS	YASARGIL
Medicon	58.56.53
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE681
AUTHORS	YASARGIL
Medicon	58.56.53

SPECIFICATIONS

[straight](#), [temporary](#)