


## YASARGIL MINI VESSEL CLIP

temporary; 3mm; opening 3.3mm

90 gms; straight



**SKU:** 244-560-03

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE681

**AUTHORS** [YASARGIL](#)

**Medicon** 58.56.53

**SPECIFICATIONS** [straight](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE681

**AUTHORS** [YASARGIL](#)

**Medicon** 58.56.53

## SPECIFICATIONS

[straight](#), [temporary](#)