


YASARGIL MINI VESSEL CLIP

temporary; 5mm; opening 3.8mm

70 gms; curved




SKU: 244-565-05

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
SPECIFICATIONS	curved , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
SPECIFICATIONS	curved , temporary