


YASARGIL MINI VESSEL CLIP

permanent; 3mm; opening 3.3mm

110 gms; straight



SKU: 244-462-03

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE700
AUTHORS	YASARGIL
Medicon	58.56.04
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE700
AUTHORS	YASARGIL
Medicon	58.56.04

SPECIFICATIONS

[permanent](#), [straight](#)