


## YASARGIL MINI VESSEL CLIP

temporary; 4mm; opening 3.6mm

80 gms; curved



**SKU:** 244-567-04

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.57.61

**SPECIFICATIONS** [curved](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.57.61

**SPECIFICATIONS** [curved](#), [temporary](#)