


YASARGIL MINI VESSEL CLIP

temporary; 4mm; opening 3.6mm

80 gms; curved



SKU: 244-567-04

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.57.61

SPECIFICATIONS [curved](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.57.61

SPECIFICATIONS [curved](#), [temporary](#)