


YASARGIL MINI VESSEL CLIP

temporary; 7mm; opening 4.4mm

70 gms; curved



SKU: 244-567-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE723

AUTHORS [YASARGIL](#)

Medicon 58.57.58

SPECIFICATIONS [curved](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE723

AUTHORS [YASARGIL](#)

Medicon 58.57.58

SPECIFICATIONS

[curved](#), [temporary](#)