


YASARGIL MINI VESSEL CLIP

permanent; 5mm; opening 4mm

110 gms; straight




SKU: 244-462-05

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE710
AUTHORS	YASARGIL
Medicon	58.56.06
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE710
AUTHORS	YASARGIL
Medicon	58.56.06

SPECIFICATIONS

[permanent](#), [straight](#)