


YASARGIL MINI VESSEL CLIP

permanent; 7mm; opening 4.6mm

110 gms; straight



SKU: 244-462-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE720
AUTHORS	YASARGIL
Medicon	58.56.07
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE720
AUTHORS	YASARGIL
Medicon	58.56.07

SPECIFICATIONS

[permanent](#), [straight](#)