


YASARGIL MINI VESSEL CLIP

permanent; 3mm; opening 3.2mm

110 gms; curved



SKU: 244-463-03

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE682
AUTHORS	YASARGIL
Medicon	58.57.03
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE682
AUTHORS	YASARGIL
Medicon	58.57.03

SPECIFICATIONS

[curved](#), [permanent](#)