


## YASARGIL MINI VESSEL CLIP

permanent; 5mm; opening 3.8mm

110 gms; curved



**SKU:** 244-465-05

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE712

**AUTHORS** [YASARGIL](#)

**Medicon** 58.57.06

**SPECIFICATIONS** [curved](#), [permanent](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE712

**AUTHORS** [YASARGIL](#)

**Medicon** 58.57.06

## SPECIFICATIONS

[curved](#), [permanent](#)