


## YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 7mm; opening 4.4mm

110 gms; curved



**SKU:** 244-467-07T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FT722T
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.43.07
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">permanent</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FT722T
<b>MATERIAL</b>	<a href="#">Titanium</a>

## AUTHORS

[YASARGIL](#)

## Medicon

58.43.07

## SPECIFICATIONS

[curved](#), [permanent](#)