


YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 7mm; opening 4.4mm

110 gms; curved




SKU: 244-467-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT722T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.07
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT722T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.43.07

SPECIFICATIONS

[curved](#), [permanent](#)