


YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 4mm; opening 3.6mm

110 gms; curved



SKU: 244-469-04T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| Aesculap | FT694T |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.43.15 |
| SPECIFICATIONS | curved , permanent |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| Aesculap | FT694T |
| MATERIAL | Titanium |

AUTHORS

[YASARGIL](#)

Medicon

58.43.15

SPECIFICATIONS

[curved](#), [permanent](#)