


YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 4mm; opening 3.5mm

110 gms; curved



SKU: 244-471-04T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT714T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.16
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT714T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.43.16

SPECIFICATIONS

[curved](#), [permanent](#)