


YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 7mm; opening 6mm

130 gms; straight



SKU: 244-475-06T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT717T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.41
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT717T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.43.41

SPECIFICATIONS

[permanent](#), [straight](#)