


YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 5mm; opening 4mm

110 gms; straight



SKU: 244-460-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT690T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.42.05
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT690T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.42.05

SPECIFICATIONS

[permanent](#), [straight](#)