


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 3mm; opening 3.3mm

90 gms; straight



SKU: 244-562-03T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.42.63
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.42.63

SPECIFICATIONS

[straight](#), [temporary](#)