


## YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 3mm; opening 3.3mm

90 gms; straight



**SKU:** 244-562-03T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                 |
| <b>TYPE</b>           | <a href="#">Vessel Clips</a>                         |
| <b>MATERIAL</b>       | <a href="#">Titanium</a>                             |
| <b>AUTHORS</b>        | <a href="#">YASARGIL</a>                             |
| <b>Medicon</b>        | 58.42.63   |
| <b>SPECIFICATIONS</b> | <a href="#">straight</a> , <a href="#">temporary</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                 |                              |
|-----------------|------------------------------|
| <b>UNIT</b>     | <a href="#">pcs.</a>         |
| <b>TYPE</b>     | <a href="#">Vessel Clips</a> |
| <b>MATERIAL</b> | <a href="#">Titanium</a>     |
| <b>AUTHORS</b>  | <a href="#">YASARGIL</a>     |
| <b>Medicon</b>  | 58.42.63                     |

## SPECIFICATIONS

[straight](#), [temporary](#)