


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 7mm; opening 4.4mm

70 gms; curved



SKU: 244-567-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.58
SPECIFICATIONS	curved , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.58

SPECIFICATIONS

[curved](#), [temporary](#)