


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 4mm; opening 3.5mm

80 gms; curved



SKU: 244-571-04T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.67
SPECIFICATIONS	curved , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.67

SPECIFICATIONS

[curved](#), [temporary](#)