


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 5mm; opening 4mm

70 gms; curved



SKU: 244-571-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.68
SPECIFICATIONS	curved , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.68

SPECIFICATIONS

[curved](#), [temporary](#)