


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 5mm; opening 3.5mm

70 gms; straight



SKU: 244-573-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	straight , temporary