


## YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 4mm; opening 7mm

70 gms; angled



**SKU:** 244-576-04T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FT227T

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.43.71

**SPECIFICATIONS** [angled](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FT227T

**MATERIAL** [Titanium](#)

## AUTHORS

[YASARGIL](#)

## Medicon

58.43.71

## SPECIFICATIONS

[angled](#), [temporary](#)