


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 4mm; opening 7mm

70 gms; angled



SKU: 244-576-04T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT227T

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.43.71

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT227T

MATERIAL [Titanium](#)

AUTHORS

[YASARGIL](#)

Medicon

58.43.71

SPECIFICATIONS

[angled](#), [temporary](#)