


YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 5mm; opening 4mm

110 gms; straight



SKU: 244-462-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT710T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.42.06
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT710T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.42.06

SPECIFICATIONS

[permanent](#), [straight](#)