


## YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 7mm; opening 4.6mm

110 gms; straight



**SKU:** 244-462-07T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FT720T
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.42.07
<b>SPECIFICATIONS</b>	<a href="#">permanent</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FT720T
<b>MATERIAL</b>	<a href="#">Titanium</a>

## AUTHORS

[YASARGIL](#)

## Medicon

58.42.07

## SPECIFICATIONS

[permanent](#), [straight](#)