


YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 4mm; opening 3.8mm

110 gms; curved



SKU: 244-463-04T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT692T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.05
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT692T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.43.05

SPECIFICATIONS

[curved](#), [permanent](#)