


YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 4mm; opening 3.6mm

110 gms; curved



SKU: 244-465-04T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.08
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.08

SPECIFICATIONS

[curved](#), [permanent](#)