


YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 3mm; opening 3.3mm

110 gms; straight



SKU: 244-460-03T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.42.03
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.42.03

SPECIFICATIONS

[permanent](#), [straight](#)