


YASARGIL VESSEL CLIP; STANDARD

temporary; 7mm; opening 7.9mm

90 gms; angled



SKU: 244-510-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.59.06

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.59.06

SPECIFICATIONS [angled](#), [temporary](#)