


## YASARGIL VESSEL CLIP; STANDARD

temporary; 12mm; opening 10.2mm

90 gms; angled



**SKU:** 244-510-12

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)


**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.59.12

**SPECIFICATIONS** [angled](#), [temporary](#)

**Schließkraft**  
**Closing Force**



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.59.12

**SPECIFICATIONS** [angled](#), [temporary](#)