


YASARGIL VESSEL CLIP; STANDARD

temporary; 5mm; opening 6.2mm

110 gms; angled



SKU: 244-511-05

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.59.85

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.59.85

SPECIFICATIONS [angled](#), [temporary](#)