


YASARGIL VESSEL CLIP; STANDARD

permanent; 15mm; opening 8.7mm

200 gms; curved



SKU: 244-401-15

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| AUTHORS | YASARGIL |
| SPECIFICATIONS | curved , permanent |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| AUTHORS | YASARGIL |
| SPECIFICATIONS | curved , permanent |