


YASARGIL VESSEL CLIP; STANDARD

temporary; 8mm; opening 5.4mm

130 gms; curved



SKU: 244-515-08

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.63.08

SPECIFICATIONS [curved](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.63.08

SPECIFICATIONS [curved](#), [temporary](#)