


YASARGIL VESSEL CLIP; STANDARD

temporary; 7mm; opening 5.7mm

90 gms; straight



SKU: 244-521-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.63.37

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.63.37

SPECIFICATIONS [straight](#), [temporary](#)