


## YASARGIL VESSEL CLIP; STANDARD

temporary; 11mm; opening 8mm

90 gms; straight



**SKU:** 244-521-11

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.63.41

**SPECIFICATIONS** [straight](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.63.41

**SPECIFICATIONS** [straight](#), [temporary](#)