


YASARGIL VESSEL CLIP; STANDARD

temporary; 9mm; opening 6.2mm

90 gms; straight



SKU: 244-523-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.63.49

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.63.49

SPECIFICATIONS [straight](#), [temporary](#)