


YASARGIL VESSEL CLIP; STANDARD

temporary; 7mm; opening 7.4mm

110 gms; straight; fenestrated Ø 3.5mm



SKU: 244-530-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.07

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.07

SPECIFICATIONS [straight](#), [temporary](#)