


## YASARGIL VESSEL CLIP; STANDARD

temporary; 12mm; opening 9.6mm

90 gms; straight; fenestrated  $\varnothing$  3.5mm



**SKU:** 244-530-12

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN


**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.58.12

**SPECIFICATIONS** [straight](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.58.12

**SPECIFICATIONS** [straight](#), [temporary](#)