


YASARGIL VESSEL CLIP; STANDARD

temporary; 9mm; opening 9.1mm

100 gms; straight; fenestrated Ø 5mm



SKU: 244-536-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.39

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.39

SPECIFICATIONS [straight](#), [temporary](#)