


## YASARGIL VESSEL CLIP; STANDARD

temporary; 5mm; opening 6.5mm

110 gms; angled; fenestrated Ø 5mm



**SKU:** 244-537-05

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN


**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.58.45

**SPECIFICATIONS** [angled](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.58.45

**SPECIFICATIONS** [angled](#), [temporary](#)