


YASARGIL VESSEL CLIP; STANDARD

temporary; 11mm; opening 5.5mm

90 gms; angled; fenestrated Ø 5mm



SKU: 244-539-11

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.61

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.61

SPECIFICATIONS [angled](#), [temporary](#)