


## YASARGIL VESSEL CLIP; STANDARD

temporary; 11mm; opening 10mm

90 gms; angled; fenestrated Ø 5mm



**SKU:** 244-540-11

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.58.71

**SPECIFICATIONS** [angled](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.58.71

**SPECIFICATIONS** [angled](#), [temporary](#)