


## YASARGIL VESSEL CLIP; STANDARD

temporary; 10mm; opening 5.8mm

110 gms; angled; fenestrated Ø 5mm



**SKU:** 244-543-10

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.41.20

**SPECIFICATIONS** [angled](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.41.20

**SPECIFICATIONS** [angled](#), [temporary](#)