


YASARGIL VESSEL CLIP; STANDARD

temporary; 10mm; opening 8.8mm

110 gms; straight; fenestrated Ø 5mm



SKU: 244-545-10

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.41.30

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.41.30

SPECIFICATIONS [straight](#), [temporary](#)