


YASARGIL VESSEL CLIP; STANDARD

permanent; 15mm; opening 7.4mm

200 gms; curved



SKU: 244-403-15

Categories: [Neurosurgery](#), [Open Surgery](#)


ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

SPECIFICATIONS [curved](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

SPECIFICATIONS [curved](#), [permanent](#)