


YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 8.7mm

200 gms; angled



SKU: 244-410-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE758

AUTHORS [YASARGIL](#)

Medicon 58.55.45

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE758

AUTHORS [YASARGIL](#)

Medicon 58.55.45

SPECIFICATIONS

[angled](#), [permanent](#)